1- Indications

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1- Indications

1A200 **GERY** foot is designed exclusively for lower limb prosthesis. It is appropriate for low activity patients (L1) weighing (carried load included) less than 100 kg for sizes 22 to 27 cm, and less than 125 kg for sizes 28 and 29 cm. During the manufacture of the prosthesis, this foot, whatever its size, can be adapted to shoes with heel height varying from 0 mm to 20 mm.

! CAUTION : 1A200 GERY foot is not appropriate for bilateral amputees.

2- Description and functions

1A200 **GERY** foot includes two main components which cannot be dismantled : the core and the cosmetic cover.

It also includes a malleolar clip used to fix the cosmetic foam.

The cosmetic cover, made of a material that does not react to moisture and UV, upholsters the structure and ensures a good position in a shoe of the convenient size. These components reduce the stresses in the residual limb during walk, allowing the foot to be flexible and comfortable.

The hole at the lower face of the foot is used as a guide mark if it is necessary to separate the big toe from the others. In this case, cut the foot up to the hole with a brand saw and widen the cut as needed.

3- Assembling and alignment

Use compatible links to assemble the foot as usual and check that the load line falls between 0 and 10 mm before the pyramid axis.

4- Dynamic adjustment

According to your patient's comments and what you notice, make the adjustments that will ensure the least tiring and most comfortable walk for the patient.

5- Installation of the malleolar clip

Cut the cosmetic foam to the desired length, position it on the foot and draw the malleolus contour.

Inside the foam mill a space that will receive the selected connector, and at the same time take care that thickness be sufficient to allow for sticking.

Clean contact surfaces and glue the malleolar clip on the foam with XC090 neoprene adhesive.

After drying, finish the outer shape as usually.

6- Maintenance



Zero maintenance.

Clean the cosmetic cover with a slightly damp cloth or sponge.

7- Recycling

If possible, remove the core from the cosmetic cover, then separate the plastic part of the core from the metal part (pyramid and sphere).

The cosmetic cover and the metal part of the core are specific wastes. The metal part is made of stainless steel. These parts must be recycled in compliance with the laws in force.

! SAFETY INSTRUCTIONS FOR THE PATIENT

- Do not let foreign bodies inside the foot.
- If the foot has been dipped or sprinkled, empty the remaining liquid out and dry the foot in the open. In case of metallic component corrosion or doubts about the behaviour of the prosthesis, have your prosthesis controlled.
- Clean the cosmetic cover with a slightly damp cloth or sponge.
- Use temperature : - 20° / + 60°C
- An inappropriate use of the foot (walking intensively, sport, frequent overload, ...) can damage components and be dangerous for you. If you discover an abnormal behaviour of your foot or changes in your foot features, contact your prosthetists.
- Changing shoes (with a different heel height) can cause discomfort and pains in the stump. It can also lead to changes in the feeling of balance.
- While going upstairs, take care to put your foot flat on the step. Do not apply your weight on the forefoot only.

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