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Fig. 14


Fig. 15

ASSEMBLY AND ADJUSTMENT INSTRUCTIONS FOR 1P50 - R

CONNECTION WITH OTHER PROSTHETIC MODULES :

Connection in upper part with any Proteor standard module for knee linking, and in lower part with any SACH or DYNASTEP foot from PROTEOR range.

USE :


For patients up to 100 kg, and for any activity level.

I : PACK CONTENTS

(See figure 1)

- A Carbon frame
- B Hydraulic unit
- C Rod
- D Upper joint
- E Lower joint
- F Hexagonal hollow screw with cylindrical head (x 4)
- G Spacers (x 4)
- G1 Plain washer
- G2 Spring Grower
- H1 Hexagonal hollow screw with cylindrical head M10x30
- H2 Hexagonal hollow screw with cylindrical head M10x45
- H3 Hexagonal hollow screw with cylindrical head M10x60
- Z tool for locking in rotation

Fig. 1

CAUTION :

- 1) When cutting and assembling, never modify the position of the green knob of the hydraulic unit.
- 2) The cutting block ref. 1 P 5070 must **IMPERATIVELY** be used for adjusting the length of the frame and that of the rod. **PROTEOR cannot be held liable if the Hydracadence frame is incorrectly cut.**

II : CHECK THE SCREW FITTINGS FOR THE FOOT

(See figures 2, 3 and 4)

Feet references 1A21 and 1A118 :

- Spring Grower **G2**
- Screw **H2** (M10x45)

Fig. 2



Fig. 3

Feet reference 1A101 :


- Spring grower **G2**
- Screw **H1** (M10x30)

Dynastep Foot 1A101 when fitted with the frame been sawn to the length relevant to feet 1A21 or 1A118 :

- Spring grower **G2**
- Screw **H3** (M10x60)
- Wooden group 1D70/1D71

Fig. 4

III : LENGTH ADJUSTMENT OF THE FRAME AND THE ROD (See figure 5)

CAUTION : Check when fitting that the piston *M* is in the neutral position. To make sure, press the end of the rod *C* with your thumb to move back the piston *M*, then release the pressure to allow the piston *M* to return to the neutral position.

1 - Fit the prosthesis to measure the dimension X (See figure 5)


X = knee joint / heel

As indicated in figure 5, fit the following parts :

- A** Complete frame with hydraulic unit
- C** Rod (not connected with the eyelet **P**)
- E** Lower joint
- J** Foot

- Be care that the lower part of the frame be stopped by the lower joint (see the contact point indicated by * on figure 5).
- Fit the tool for locking in rotation **Z**.
Z2 threaded shank will get through the **C** rod head axis hole.
- Place under the heel a wedge **L** corresponding to the height of the shoes instep usually worn by the patient.
- Keep the assembly vertical and measure the distance **X** as indicated in figure 5.

Fig. 5

2 - Calculating length Y to be cut (See figure 5)

. Measure on the patient his knee/ground distance. ***This measurement must be made without shoes. For the knee, the measure will be considered as beginning at the physiological joint axis.***
KNEE / GROUND MEASURED ON A PATIENT = W

- Calculate the cut length Y

This value is the result of the following formula :


$$Y = X - W + 4\text{mm}$$

3 - Frame and rod cutting, connection holes

***CAUTION! : Do not adjust the green knob of the hydraulic unit.
Keep in place the tool for locking in rotation.***

. Fit the frame A / hydraulic unit B / rod C in the cutting block as shown in figure 6, taking care to adjust the

Fig. 6

- Use a hacksaw with a blade for cutting steel to cut the rod C (see figure 6) along surface S of the cutting block (see fig. 6), then cut the frame A along surface T of the cutting block.
- With a $\varnothing 7$ mm drill, bore 4 connection holes using the guide holes R of the cutting block (see figure 6).
- Remove the { frame / hydraulic unit / rod } assembly from the cutting block and burrs.

IV : ASSEMBLING (See figures 7, 8, 9 and 10)

1 - Fitting the ankle joint

Detach the eyelet from the joint assembly (See figure 7)

- Undo the screw N
- Drive out the pin O
- Remove the eyelet P


Fig. 7

Fit the eyelet to the rear rod (See figure 8)

CAUTION! : Keep in place the tool for locking in rotation.

- Screw
- Do not tighten the lock nut V.



Fig. 8

Fit the ankle joint onto the frame (See figure 9)

Fig. 9


- 1) Slide the joint **E** into the frame **A**
- 2) Position the spacers **G** in the frame holes
- 3) Tighten the screws **F** after having coated them with average strength adhesive ref. XC047.
Advised tightening torque value : 11 Nm

Fit the rod and the ankle joint together (See figure10)

- 1) Fit the rod **C** + eyelet **P** assembly into the rear shackle of the joint **E** (see figure10).
- 2) Refit the pin **O** making sure that one of the holes * is lined up with the tapped hole for the screw **N** (see figure 10).
- 3) Tighten the screw **N** after having coated it with average strength adhesive PROTEOR ref. XC047.

Fig. 10


CAUTION :
Do not tighten the eyelet lock nut yet.

2 - Feet assembling (See figure 11)

CAUTION : Do not tighten the eyelet lock nut yet. Keep in place tool Z.

Fig. 11

Advised tightening torque
value : 22 Nm

V : CHECKING (See fig.12)

a) After fitting the socket, check that when flexed, this socket comes into contact with the hoop **A1** of the frame before the hydraulic unit **B** comes up against the inside **A2** of the frame.

b/ Z1 tapped heads will be supported by frame A edges.

c) Tighten the lock nut **V**.

d) Dismantel the tool for locking in rotation **Z**.

VI - ADJUSTMENT (See figures 13, 14 and 15)

Getting the perpendicularity of the hydracadence unit.

a) Keep the hydracadence unit vertical and make a mark with a pencil at 12 mm in front of the articular center of the ankle (figure 13).

b) Use the heel height adjustment knob **B4** to modify the place of the articular center of the knee (figure 14).

The articular center of the knee must be on the vertical line crossing at the mark situated 12 mm in front of the ankle joint.

c) The prosthesis being still on an horizontal area, check with a level that the upper face is perfectly horizontal. If necessary, modify the adjustment with the knee piston **B1** (figure 15).

To do so :

- Loosen the nut **B2**.
- Turn the piston **B1** with the special wrenches reference 5F01 and 5F04.
- When adjusted as required, retighten the nut **B2**.

IMPORTANT

After this adjustment, remove the smear on the small holes located at the upper level of the piston rod **B1**.

d) Pace adjustment

To adjust the pace, use the red based knob **B3**.

- Turn clockwise to lock completely, then give 3 and half wrench turns or 7 half-turns to open. Never exceed this limit because it could damage the rubber gasket.

This is the limits of the adjustment. This adjustment will decrease or increase the entry of the fluid, creating resistance or not.

e) Heel height adjustment

To change the height of the heel without altering the perpendicular setting of the prosthesis (e.g. change of foot wear), use the green knob **B4**.


Fig. 13

Fig. 14



Fig. 15

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