
NOTICE D'ASSEMBLAGE ET DE REGLAGE POUR GENOU 1P340
2/17 —> 5/17
-
INSTRUCTIONS FOR ASSEMBLY AND ADJUSTMENT OF KNEE 1P340
6/17 —> 9/17
-
AUFBAU- UND EINSTELLANLEITUNG FÜR KNIE 1P340
10/17 —> 13/17
-
INSTRUCCIONES DE MONTAJE Y DE AJUSTE PARA RODILLA 1P340

NABTESCO CORPORATION
35, Uozakihama-machi-Higashinada-ku
Kobe 658-0024 JAPAN
PHONE : +33 (0)3 80 78 42 10 - FAX : +33 (0)3 80 78 42

NABMIC B.V
Brouwersrtaat 34 - 2984 AR Ridderkerk
THE NETHERLANDS
PHONE : 0180-530-590 - FAX : 0180.530.591

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INSTRUCTIONS FOR ASSEMBLY AND ADJUSTMENT FOR 1P340 KNEE

I - ALIGNMENT SETTING

When assembling the trans-femoral prosthesis with 1P340 knee, set the basic alignment as follows (fig. 1a and 1b) :

a) Alignment setting in frontal plane

Set the alignment so that the load line of the body weight goes through the centre of the knee and the centre of the heel (fig. 1a).

b) Alignment setting in sagittal plane

Set the alignment so that the load line of the body weight goes 5 mm in front of the centre of the upper front bar and through the mid-point between heel and toe-break (fig. 1b).

c) Additional settings

After these basic alignment settings, we recommend :

- For very active patients : further backwards adjustment of the alignment : between 0 and 5 mm in front of the centre of the upper front bar
- For less active patients : further forwards adjustment of the alignment : between 5 to 10 mm in front of the centre of the upper front bar.

Note : During this adjustment, pay attention to the following :

- *Do not set the alignment beyond the above figures. Otherwise this will create unreasonable forces against the knee and may cause trouble or damage to the knee.*
- *Set the alignment when the knee joint is fully extended.*


Fig. 1a
Basic alignment :
Frontal plane


Fig. 1b
Basic alignment :
Sagittal plane

II - LOOSENESS ADJUSTMENT AT KNEE AXIS

Some looseness may appear in the joint after a few months of use. It can be corrected by turning toothed knobs **A**, then **B**, then **C** (fig. 2).


Fig. 2

a) To measure the tightening amount of the toothed knobs, make a mark on the root thread of knob **A** which is in contact with the clamping screw (Fig. 3).


Fig. 3

b) Remove the clamping screw using a 1,5 mm hex wrench (Fig. 4).


Fig. 4

c) With a special tool used to tighten toothed knobs, tighten Knob **A** by a single thread (fig. 5), then control the degree of looseness. If some looseness remains, apply further tightening.

Caution: Over-tightening will increase resistance to knee flexion and extension (Fig. 5).


Fig. 5

d) If the above indicated tightening of knob **A** solves the looseness problem, then it is not necessary to adjust knobs **B** and **C**.

If some looseness still remains, use the same procedure to adjust knobs **B**, then **C**.

e) Make knee flexions and extensions to check flexibility, then install the clamping screws (Fig. 6). If the resistance to flexion or extension is too high, loosen toothed knobs **C**, then **B**, then **A**.


Fig. 6

Apply Loctite adhesive on each locking screw.

III - BATTERY REPLACEMENT

1. Remove the battery connector off the circuit board (Fig. 7).
2. Loosen the fastening screws off the outer clamp, then pull out the pipe.
3. Remove the outer clamp off the frame, then the fastening screw, then pull out the inner clamp together with the battery case (Fig. 8).
4. Using a coin, turn the cap of the battery case, then remove the battery (Fig. 9).
5. Place a new battery in the battery case. The battery pipe must be in the lit of the inner clamp. Otherwise the battery will not go fully into depth. (Fig. 10).
6. Set the cap of the battery case, then fasten it.
7. Set the inner clamp in the frame. Check that the battery connector is on the front face of the circuit board (Fig. 8).
8. Install the fastening screw and tighten it so that it does not protrude the frame surface.
Caution : Over-tightening may cause damage to the battery case and perhaps to the battery itself.
9. Set the outer clamp.
10. Plug the battery connector in the circuit board. Using tweezers can ease this operation.


Fig. 7 : Removing the battery connector


Fig. 8 : Removing the clamp


Fig. 9 : Removing the battery cap

Fig. 10 : Removing the battery

IV - REPLACEMENT OF RUBBER STOPPERS

Rubber stoppers are located at two different places : one for extension and the other for flexion.
They will be more or less worn out according to the frequency of use.

However, we recommend to replace them each year.

The replacement must be done as follows :

a) Extension rubber stopper (Fig. 11)

- The knee being bent, use a minus driver to pick up and take out the rubber stopper.
- Set a new rubber stopper.

Caution : Even if the installation of this stopper was hard, check that it has been made correctly.

b) Flexion rubber stopper (Fig. 12)

- The knee being completely extended, use a minus driver to pick up and take out the rubber stopper.
- Set a new rubber stopper.

Caution : Even if the installation of this stopper was hard, check that it has been made correctly.


Fig. 11 :
Extension
stopper


Fig. 12 :
Flexion
stopper

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