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

### 1 – DESCRIPTION AND FUNCTIONS

1A101 DYNASTEP foot includes a polyurethan foam cosmetic structure **1** moulded on two composite carbon blades **2** and **3**. The original shape of the blades ensures great walking comfort and a slight energy return that allows to spare energy during walking and to increase the patient's possible walking distance. Blade **2** is split to help inversion-eversion motions. Big toe is apart.


Thanks to DYNASTEP, trans-tibial and trans-femoral amputees will enjoy a comfort matching to their activities and their physical abilities in their everyday life.

DYNASTEP comes with its malleolar clip to help the set up of the cosmetic shell.

### 2 – CONNECTION

DYNASTEP foot can be used with any ankle connector **A** that includes a M10 screw (**A1**).

### 3 – SPECIFICATIONS

DYNASTEP is designed exclusively for lower limb prosthesis. It is appropriate for moderate activity patients weighing (carried load included) less than 80 kg for sizes 22/23/24/25 and less than 100 kg for sizes 26/27/28/29 cm.

### 4 – ALIGNMENTS

#### A) Heel height

Due to its foot arch, DYNASTEP is appropriate for shoes with 10 mm heel height. However it can be used with shoes whose heel height ranges from 5 to 15 mm if it is combined with an adaptor providing equinus adjustment.


*Note: Shoes with heel height above 10 mm tend to increase energy return and shoes with heel height under 10 mm tend to decrease energy return.*

#### B) Load line

Align the foot so that the load line goes through the posterior third of the foot.

Notes:

- Increasing equinus provides a more dynamic alignment. It increases energy return at step completion.
- If the foot seems too long during transfer to the forefoot or if heel support is too flexible, you have to position the foot behind the socket.
- If heel support seems too hard or if the foot seems too short during transfer to the forefoot, you have to position the foot ahead of the socket.

### 5 – INSTALLATION OF THE MALLEOLAR CLIP

Cut the cosmetic foam to the desired length, position it on the foot and draw the malleolus contour.

Inside the foam mill a space that will receive the selected connector, and at the same time take care that thickness be sufficient to allow for sticking.

Clean contact surfaces and glue the malleolar dip to the foam with XC090 neoprene adhesive.

After drying, finish the outer shape as usual.

### 6 – RECYCLING

The DYNASTEP foot includes a polyurethan foam cosmetic structure moulded on two composite blades. It must be recycled according to the laws in force.

Part for the patient

### ADVICE FOR USE, MAINTENANCE AND SAFETY

Clean the cosmetic cover with a slightly damp cloth or sponge.

Use temperature : - 20°C / + 60°C.

Do not let foreign bodies or dust inside the foot.

If the foot has been dipped or sprinkled, empty the remaining liquid and dry the foot rapidly.

A use of the foot that does not adhere to your prosthetist's recommendations can damage foot components.

In case of an abnormal behaviour of your foot, or changes in your foot features, or if the foot sustained a hard shock, contact your prosthetist.

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